




# Applanation Tonometer AT-1

Aiming at new levels in quality



# Applanation Tonometer AT-1

## Applicable tonometer base

For SM-70N



S12-13

For SM-30N and SM-10N



S12-14

For SM-4N and SM-2N



S12-15

For 300XL



S12-16



## Specifications

Model	AT-1				
Model of slit lamp	SM-90 SM-90N	SM-70 SM-70N	SM-30,SM-10 SM-30N,SM-10N	SM-4,SM-2 SM-4N,SM-2N	300XL
Tonometer main unit	S12-12				S12-12 XL
Tonometer base	—	S12-13	S12-14	S12-15	S12-16
Type of AT-1	AT-1	AT-1 Type 7	AT-1 Type 3	AT-1 Type 4	AT-1 Type XL
Measurement method	Depression type (observed with left eye on the slit lamp)				
Depression area	7.354mm <sup>2</sup>				
Diameter of measurement depression area	3.06mm				
Measurement range	0 ~ 80mmHg				
Minimal scale	2mmHg				
Measurement error	±1mmHg				
Equation for conversion from mmHg to kpa	1mmHg = 101325/760Pa ≈ 133.3Pa 100Pa = 0.1kPa				

● Design and specifications are subject to change as improvements are made to the product.

**TAKAGI SEIKO CO.,LTD.**

330-2 IWAFUNE, NAKANO-SHI, NAGANO-KEN, 383-8585, JAPAN  
TEL.+81-269-22-4512 FAX.+81-269-26-6321  
URL:http://www.takagi-j.com E-mail:info@takagi-j.com

**CE**